


KOBOT™

Nitinol Stone Retrieval Tipless Basket


- In line with the concept of Cobot, a new range of Stone Retrieval Baskets. These want to be a tool diversified, effective and integral to the work of the surgeon and his assistants and to maximize performance.
- A innovative and ergonomic Criterium Handle with assisted closing system, helps in catching, retain and remove stone fragments.



KOBOT™


Nitinol Stone Retrieval Tipless Basket

ROCAMED 





KOBOT TIPLESS

- 4 wire basket





➤ CRITERIUM HANDLE SYSTEM

Available with assisted closing System



KOBOT SIX

- 6 wire basket




KOBOT FILTER

- 16 wire filter

➤ FLEXI TIP IMPROVEMENT

For a better scope deflection



KOBOT FIBER

- 4 wire basket with filter and coaxial lumen for laser fiber (200 μm)

All supplied with Tuohy Borst Connector





Ordering Informations


Product code	Size Fr	Wires	Working Length cm	Basket Length mm	Description
ROSZ1930ST	1.9	4	120	12	Kobot Tipless with criterium handle
ROSZ1940ST	1.9	4	120	12	Kobot Tipless - Criterium handle with assisted closing system
ROKT1800ST*	1.8	16	120	14	Kobot Filter with criterium handle
ROKB3000ST*	3.0	16	120	14	Kobot Fiber (200μm) with criterium handle
ROKS1700ST*	1.7	6	120	10	Kobot Six (200μm) with criterium handle

* available soon

For more information visit www.rocamed-urology.com

 9 avenue Albert II
MC 98000 Monaco

 (+377) 97 98 42 43
(+377) 92 05 61 50 (fax)

 customer@rocamed.eu

 www.rocamed-urology.com