

Key benefits

- » Gently unfolds the urethra
- » Excellent and long-lasting lubrication properties
- » Protects against injuries
- » Antiseptic action*
- » Local anaesthetic action**
- » High viscosity, no flow back from urethra***
- » Free of parabens, PVC and latex
- » Single hand use




* Cathejell products with Chlorhexidine
 ** Cathejell products with Lidocaine
 *** Montavit, internal data on file


Application of Cathejell

All Cathejell versions are perfectly suited for the operating theatre and the doctor's office

Attach the applicator tip



Instillation into the urethra



601.091/042018



SIMPLY CLEVER SINGLE HAND USE

Cathejell

Sterile & clear lubricant gel



Cathejell Mono

Cathejell C


Cathejell Lidocain C

Cathejell Lidocain



Cathejell products

- Cathejell Lidocain C**
Lubricant + anaesthetic + antiseptic
Lidocaine and Chlorhexidine
Medical device class III
- Cathejell Lidocain**
Lubricant + anaesthetic
Lidocaine
Medical device class III
- Cathejell C**
Lubricant + antiseptic
Chlorhexidine
Medical device class III
- Cathejell Mono**
Lubricant
Plain lubricant gel
Medical device class IIa



Cathejell design
Single hand use and soft instillation

First steam sterilised
single use catheter lubricant worldwide since 1986

Production in Austria
in compliance with ISO 13485



All Cathejell versions

- » Are available in 12.5 g and 8.5 g syringes for single use
- » Can be prolonged by sterile applicator tip (Cannula)

	Cathejell Lidocain C	Cathejell Lidocain	Cathejell C	Cathejell Mono
Contents (1 g gel contains)	Sterile, water-soluble, clear lubricant gel	Sterile, water-soluble, clear lubricant gel	Sterile, water-soluble, clear lubricant gel	Sterile, water-soluble, clear lubricant gel
Indications	20 mg Lidocaine Hydrochloride, 0.5 mg Chlorhexidine Dihydrochloride	20 mg Lidocaine Hydrochloride	0.5 mg Chlorhexidine Dihydrochlorid	Plain sterile gel
Indications	Lubricant for » Catheters » Endoscopes » Probes » Tubes » Tracheal intubation » Other medical instruments	Lubricant for » Catheters » Endoscopes » Probes » Tubes » Tracheal intubation » Other medical instruments	Urethral lubricant for » Catheters » Endoscopes » Probes » Tubes » Other medical instruments	Lubricant for » Catheters » Endoscopes » Probes » Tubes » Tracheal intubation » Other medical instruments
Anaesthetic	+	+		
Antiseptic	+		+	
Clear for instruments	+	+	+	+
Lubricating	+	+	+	+
Prevention of infections	+		+	
Sterilisation procedure	121° C steam sterilised			
Packages	1 x 12.5 g 5 x 12.5 g 25 x 12.5 g 1 x 8.5 g * 5 x 8.5 g 25 x 8.5 g *available only for Cathejell Lidocain and Cathejell C			
Available accessory	Cathejell Cannula - Sterile applicator tip			
	Medical device class III	Medical device class III	Medical device class III	Medical device class IIa
	Adults, adolescents and children over 2 years	Adults, adolescents and children over 2 years	Doctor should be consulted for children under 5 years	No age restrictions
	Do not use Cathejell in the eyes, middle ear or wounds.			