

Order online  
minzehealth.com



## HOMEFLOW

AUTOMATED UROFLOW VOIDING DIARY AT HOME

### MINZE HOMEFLOW

#### AUTOMATED UROFLOW VOIDING DIARY AT HOME



##### On any toilet

Use on any normal home toilet



##### Multiple flows

Multiple flows, collected at home



##### Voiding diary app

Integration of uroflowmetry & voiding diary in one test







##### Automated processing

Automated processing & visualization of voiding diary




www.minzehealth.com  
info@minzehealth.com

Follow us on:  
  



## HOMEFLOW - UROFLOWMETRY ONLY

Homeflow allows uroflowmetry from home. The uroflowmeter registers volume, time and flow automatically. When the uroflowmeter registers a void, the patient is prompted in the app to complete basic questions such as type of void and urge.




## HOMEFLOW - WITH VOIDING DIARY

Homeflow can be used as a fully automated voiding diary. Volume and flow are registered automatically. The patient can add liquid intake or leakage episodes through the Voiding Diary App.

### Automatic patient identification



### Natural voiding position for men, women and children



### Processed at clinician portal

