


ROCAMED 

KOBOT™ MITT

Nitinol Stone Retrieval Tipless Basket



» » » ENDO-UROLOGY

www.rocamed.com



KOBOT™ MITT

Nitinol Stone Retrieval Tipless Basket

ROCAMED 




➤ **EASY CATCH**
Retention strength better than a tipless basket



➤ **CAPTURE, RELEASE, REPOSITION,**
and remove multiple stones and fragments in lower calix



➤ **KOBOT HANDLE TOUHY BORST**
Better grip





Ordering Informations

Product code	Size Fr	Wires	Working Length (cm)	Basket Ø (mm)	Basket Length (mm)	Description
ROKM1900ST	1.9	3	120	11	15	Kobot MITT Stone extractor

MKG-FY-MITT-EN-00
Designed by FOCALYS. Illustrations in this document are non-contractual.
ROCAMED reserves the right to change product specifications without notice.

For more information visit www.rocamed.com

 Rocamed SAM
9 avenue Albert II
98000 Monaco

 (+377) 97 98 42 43

 info@rocamed.com

 www.rocamed.com