

**SINGLE-USE DIGITAL
FLEXIBLE URETEROSCOPE**

PU3033A

7.5 Fr

Precise & Extraordinary

Reduced Outer Diameter

Irrigation Up To


19%

147%

Both Compared with PU3022A in 12/14 Fr Sheath

In 10/12 Fr sheath, irrigation is 280% compared with PU3022A


In 11/13 Fr sheath, irrigation is 179% compared with PU3022A




Optical System	
Field of View	90°
Direction of View	0°
Depth of Field	3-50 mm
Image Resolution	5 Lp/mm
Insertion Portion	
Controllable Portion	≥270° Up ≥270° Down
Nominal Width of Insertion Cord	Optics 7.5 Fr Tip 7.5 Fr
Working Length	650 mm
Working Channel	3.6 Fr



9.5 Fr ID Sheath
Adaptable



Lower
Complication Rate



Higher
Stone-Free Rate

PUSEN |

Zhuhai Pusen Medical Technology Co., Ltd.

5/F Building 1, No.33, Keji San Road

High-tech Zone, Tangjiawan Town,

519085 Zhuhai, Guangdong, People's Republic of China

www.pusenmedical.com

Rev.202104

For medical professionals only

Pusen

PUSEN Medical

@PusenMedical