PUSEN

Single Use Flexible Video Ureteroscope

With Surgeon Controlled Pressure System

PU3033H | PU3033AH

Precision Controlled Suction Innovative Urological Solutions





Integrated Suction Port

- Facilitate intrarenal pressure and temperature maintaining a safe threshold
- Optimal irrigation flow during laser activation
- Low-risk of postoperative infections and other complications



Straight Channel



Integrated T-Connector



Ergonomic Handle



Kink Resistant

Opti	ical Systen	n
Field of View		120°
Direction of View		0°
Depth of Field		3-50 mm
Inser	tion Portic	on
Controllable Portion		270° Up 270° Down
Maximum Insertion Portion Width		2.7 mm
Diameter of Insertion Portion		7.5 Fr
Distal End Diameter		7.5 Fr
Working Length		650 mm
Minimum Instrument Channel Width		3.6 Fr
Additi	onal Featu	res
Custom Buttons	PU3033I	H ○ PU3033AH ●
Irrigation Tube Included	PU30331	

PUSEN

Zhuhai Pusen Medical Technology Co., Ltd.

5/F Building 1, No.33, Keji San Road High-tech Zone, Tangjiawan Town, 519085 Zhuhai, Guangdong, People's Republic of China (86) 400-688-2023





